- (FILE 'HOME' ENTERED AT 11:47:18 ON 19 OCT 2004) FILE 'CA' ENTERED AT 11:47:35 ON 19 OCT 2004
- L1 47915 S EXPERIMENT? (3A) (DESIGN? OR SETUP OR PROCESS)
- L2 2717 S L1 AND (COMPUTER OR PROCESSOR OR MICROPROCESSOR OR AUTOMAT?)
- L3 283 S L1 AND (COMPUTER OR PROCESSOR OR MICROPROCESSOR OR AUTOMAT?) (5A) CONTROL?
- L4 54 S L1 (5A) (COMPUTER OR PROCESSOR OR MICROPROCESSOR OR AUTOMAT?) (5A) CONTROL?
- L5 114 S L3 AND (OPTIMI? OR SYNTHE? OR LIQUID OR SOLUTION OR FLUID OR SAMPLE)
- L6 2362 S L1 AND FACTORIAL
- L7 774 S L6 AND (TEMPERATURE OR L2)
- L8 346 S L7 AND (COMPOSITION OR CONCENTRATION)
- L9 44 S L8 AND (SAMPLE OR SOLUTION OF LIQUID OR FLUID)
- L10 302 S L8 NOT L9
- L11 2 S L10 AND (EXPERT? OR PHASE DIAGRAM)
- L12 39159 S (COMPUTER OR PROCESSOR OR MICROPROCESSOR OR AUTOMAT?) (5A) CONTROL?
- L13 2283 S (COMPUTER OR PROCESSOR OR MICROPROCESSOR OR AUTOMAT?) (5A (HEATER OR HEATING OR TEMPERATURE) (6A) (CONTROL? OR CHANG? OR VARIATION)
- L14 38 S L2 AND (SAMPLE OR SOLUTION OR LIQUID OR FLUID) (6A) (CHARACTERI? OR PROPERTY)
- L15 33 S L13 AND (SAMPLE OR SOLUTION OR LIQUID OR FLUID) (10A) (CHARACTER? OR PROPERTY)
- L16 2030 S L12 AND L13
- L17 251 S L16 AND (COMPOSITION OR CONCENTRATION)
- L18 83 S L17 AND (PH OR FLUORESC? OR COLOR? OR SHEAR? OR VISCOS? OR CONDUCTIVITY OR RESISTIVITY OR RESISTANCE)
- L19 84826 S (PHASE OR PROPERTY OR PARAMETER) (3A) (DIAGRAM OR GRAPH OR PLOT OR GRAPHIC OR DISPLAY?)
- L20 4 S L17 AND L19
- L21 85 S L12-13 AND L19
- L22 431 S L4-5, L9, L11, L14-15, L18, L20-21
- L23 328 S L22 NOT PY>1998
- L24 294 S L23 NOT(CEMENT OR CARS OR ALLOY OR FLASH PHOTOLYSIS OR MHD OR SOIL)
- L25 276 S L24 NOT (ATOMIC ABSORPTION OR NUCLEAR SPIN OR GAS MIXTURE OR MACERI? OR ORE OR DIAMOND OR DOT GLASS OR COAL)
- L26 262 S L25 NOT(DIFFRACT? OR COMBUST? OR EXTENSIBLE OR SEDIMENT OR MICROGRAV? OR IMAGE STICKING OR STORAGE RING)
- L27 248 S L26 NOT(ATOMIZ? OR SEEBECK OR NEON OR SAUSAGE OR DTA OR RARE EARTH OR ZONE REFINING OR FLAME)
- L28 14 S L26 NOT L27
- L29 1 S L28 AND DYE
- L30 225 S L27 NOT(SLUDGE OR PHOTONUCL? OR CRYSTALLIZER OR POWDER INJECT? OR CRYSTAL GROWTH OR FALLING DROP OR CAPSULE OR CRYOGEN? OR PHOTOCONDUCT? OR STEEL MANUFACT? OR MINERAL)
- L31 23 S L27 NOT L30
- L32 1 S L31 AND MATERIALS/TI

- 201 S L30 NOT (WATER SOFTENING OR CERAMIC OR STRIPPING OR COKE OR L33 ROTATING HEAD OR HOLOGR? OR THERMIONIC OR PZT OR SINK OR NOISE OR FUEL GAS)
- 182 S L33 NOT (CYTOM? OR KIDNEY OR RBMK OR DIELECTRIC OR NATURAL GAS OR L34 MUSCLE OR NUCLEAR REACTOR OR PYROLYSIS)
- L35 184 S L29, L32, L34
- => d bib,ab 1-184
- ANSWER 19 OF 184 CA COPYRIGHT 2004 ACS on STN L35
- AN 126:148384 CA
- Optima: a Windows-based program for computer-aided optimization of TIcontrolled-release dosage forms
- Lu, D. Robert; Abu-Izza, Khawla; Chen, Wei ΑU
- Coll. Pharmacy, Univ. Georgia, Athens, GA, 30602, USA CS
- Pharmaceutical Development and Technology (1996), 1(4), 405-414 SO
- The purpose of this work was to develop a computer program that assists AΒ optimization of controlled-release devices, both visually and math., using response surface methodol. (RSM). A Windows-based computer program, Optima, which interactively implemented a no. of subroutines for the optimization procedure, was developed. Optima is an integrated, user-friendly, and graphically oriented program for pharmaceutical dosage form optimization. Central composite design is implemented in the program. First- and second-order models contg. up to five variables can be fitted to the data. The user can also choose between linear and exponential individual desirability functions, and use them to construct an overall desirability function that combines all the response variables in a single response. The program can predict the optimum levels of exptl. variables, with respect to individual responses and/or the overall desirability. Optima has been successfully used in the development of sustained-released AZT-loaded microspheres. During the optimization process, three exptl. variables were investigated and four responses were measured. The exptl. design was a central composite design that was generated by the program. response values were used by the program to calc. the individual desirability functions, which were then combined into an overall desirability function. The individual responses as well as the overall desirability function were optimized by fitting to a second-order polynomial equation. The response surface were generated and optimum levels of the exptl. variables were predicted. The obsd. responses of the optimized formulation were very close to those predicted by Optima. The program proved to be a very useful, integrated tool for optimization of the controlled-release microspheres.
- ANSWER 20 OF 184 CA COPYRIGHT 2004 ACS on STN L35
- AN126:48820 CA
- Real-time, sensor-based computing in the laboratory TI
- Badmus, O. O.; Fisher, D. Grant; Shah, Sirish L. AU
- Univ. Alberta, Edmonton, AB, T6G 2G6, Can. CS
- Chemical Engineering Education (1996), 30(4), 280-285, 289 SO
- Specific control applications using personal computers, com. process AB

computer process-computer interfaces, simple processes, and Matlab/Simulink software are described in enough detail that they can be replicated by others. The real purpose of this paper, however, is to demonstrate the importance of real-time, sensor-based (RTSB) computing and how it can be easily and effectively integrated into university student labs. The paper describes: the exptl. processes used in our undergraduate computer process control lab., the process instrumentation and process-computer interface, the computer and communication systems, and typical software (operating system, application, and process control).

- L35 ANSWER 24 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 124:317822 CA
- TI Automation and optimization of glycine synthesis
- AU Fauduet, Henri; Nikravech, Mehrdad; Porte, Catherine
- CS Laboratoire de Chimie Bioorganique et Analytique, URA 499, BP 6759, 45067, Orleans, Fr.
- SO Process Control and Quality (1996), 8(1), 41-53
- Glycine was synthesized in an aq. medium through ammonolysis of monochloroacetic acid with hexamethylenetetramine as a catalyst. The automation of the process was performed on a lab. scale with the Logilap®system. This system allowed to optimize the influence of temp. and pH on the yield of the reaction. Research on the optimum yield was carried out using the exptl. design method. The first-order design was not suitable but the second-order design allowed to calc. the best conditions for obtaining the optimum yield from exptl. results. The correlation design and the three-dimensional space diagram led us to infer that a yield superior to 93% could be obtained when the reaction was achieved with a temp.
- L35 ANSWER 31 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 123:94695 CA
- TI Low temperature calorimetry for weakly interacting molecular systems
- AU Matsuo, Takasuke
- CS Dep. Chemistry Microcalorimetry Res. Center, Osaķa Univ., Osaka, 560, Japan
- SO Pure and Applied Chemistry (1995), 67(6), 911-18
- With development of small-sample calorimeters and their computer control, adiabatic low temp. calorimetry has become one of the diverse exptl. methods for studying more general low temp. properties of matter. A device to increase the efficiency of the calorimetric expt. is described. Calorimetric studies of deuterium-induced phase transitions, for which recent development of the exptl. technique was essential, are reviewed. Substances discussed are ammonium hexachlorometallates, rubidium hydrogen selenate, 5-bromo-9-hydroxyphenalenone, and their deuterated analogs.
- L35 ANSWER 35 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 122:197707 CA
- TI Automation of phase diagram recording
- AU Rouse, J.; Adamy, S. T.; Mehrteab, A.; Broze, G.

Technol. Cent., Colgate-Palmolive Co., Piscataway, NJ, 08855, USA CS SO

Journal of the American Oil Chemists' Society (1995), 72(1), 37-42

- An automated titrn. system was developed for generating data to AΒ construct phase diagrams, which are extremely useful in the development of personal and household products. The authors describe the system and how it can be used to perform dual titrn. A clear microemulsion sample is titrated with an oil until the dispersion turns cloudy (defined to be a transmittance < 90%). This mixt. is then dosed with a certain quantity of cosurfactant, more than enough to clear the mixt. The sample is again titrated with oil. This process continues until the sample no longer clears upon adding cosurfactant. The resulting measurements of oil uptake can be used to characterize the boundaries of the L1 or oil-in-H2O microemulsion region of the phase space. Expts. for up to sixteen samples can be performed, each having individual setup and operating instructions. Features include completely automated operation, computer-controlled 2-speed mixing, viscosity detection at the end-point condition, and the storage of results in a computerized format.
- ANSWER 38 OF 184 CA COPYRIGHT 2004 ACS on STN L35

ΑN 122:10925 CA

- Online copolymer composition control in the semicontinuous emulsion TIcopolymerization of ethyl acrylate and methyl methacrylate
- Leiza, Jose R.; de la Cal, Jose C.; Meira, Gregorio R.; Asua, Jose M. ΑU

Dep. Quim. Apl., Univ. Pais Vasco, San Sebastian, 20080, Spain CS

- Polymer Reaction Engineering (1993), Volume Date 1992-1993, 1(4), 461-SO
- A closed-loop strategy for copolymer compn. control in a semicontinuous AΒ emulsion polymn. system is presented. This strategy is based on a nonlinear adaptive plus proportional-integral controller that calcs. the flow rate of the more reactive monomer to be added into the reactor to produce a copolymer of a given compn. The nonlinear adaptive part of the controller is based on a simplified math. model of the process that includes an online adjustable parameter. The controller was checked by computer simulation and also exptl. verified during the semicontinuous emulsion copolymn. of Et acrylate and Me methacrylate carried out in a computer controlled exptl. setup using both purified and tech. grade monomers.
- L35 ANSWER 40 OF 184 CA COPYRIGHT 2004 ACS on STN

AN 121:159849

A rapid routine instrument for detecting nascent sedimentation and ΤI creaming in emulsions and suspensions

ΑU Meunier, Gerard

- Formulation, Castanet-Tolosan, F31320, Fr. CS
- Comunicaciones presentadas a las Jornadas del Comite Espanol de la SO Detergencia (1994), 25, 579-93
- TURBISCAN is equipped with a transmission/reflexion IR light detector AB (860 nm) which scans the entire height of a flat-bottomed tube. measures a large scale of turbidity (from abs. clearness to max. opaqueness) throughout the sample. It can perform 2000 acquisitions

per scan in < 3 s (i.e. 1 acquisition each 40  $\mu$ m). TURBISCAN is monitored by a personal computer. It can **automatically control** up to 100 pre-programmed readings, processed through WINDOWS (smoothing, superposition, zoom, shifting, etc.). Used for small vols. (7 mL) of concd. and opaque emulsions or suspensions, it detects nascent demixing phenomena (such as sedimentation, creaming or phase sepn.) over 4 times quicker than the naked eye, with better accuracy and reproducibility. Moreover, TURBISCAN is non-destructive, non-intrusive and can be fully automated. TURBISCAN is a high-performance tool for studying stability, soly., and **phase diagrams**.

- L35 ANSWER 44 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 120:109727 CA
- TI Correlations for viscosity of kraft black liquors at low solids concentrations
- AU Zaman, Abbas A.; Fricke, Arthur L.
- CS Dep. Chem. Eng., Univ. Florida, Gainesville, FL, 32611, USA
- SO AIChE Journal (1994), 40(1), 187-92
- The kinematic viscosities of several kraft black liquors from a 2-level, 4-variable, factorial-designed expt. for pulping slash pine were detd. for solids concns. from 10 to 50% and temps. ≤ 80° by glass capillary methods. The 4 pulping variables were cooking time, cooking temp., sulfidity, and effective alkali. Relationships between temp. and kinematic viscosity were developed by using free vol. and abs. rate theories. The results from these 2 methods were compared and discussed. A reduced variables method for dil. polymer solns. was used to correlate the viscosity with the combined effects of temp. and solids concn. The purpose of the study was to evaluate the utility of various fundamentally based models for correlating viscosity data of black liquors as a function of temp. and concn. of nonvolatile components in the region in which the liquors behave as Newtonian fluids.
- L35 ANSWER 50 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 119:170326 CA
- TI Computer-controlled system designed to measure photodegradation of photochromic compounds
- AU Dubest, R.; Levoir, P.; Meyer, J. J.; Aubard, J.; Baillet, G.; Giusti, G.; Guglielmetti, R.
- CS Inst. Topol. Dyn. Syst., Univ. Paris, Paris, 75005, Fr.
- SO Review of Scientific Instruments (1993), 64(7), 1803-8
- The authors describe a new, entirely computer-controlled app. designed for degrdn. studies of photochromic compds. Photodegrdn. can be performed in three different ways using, in addn. to the usual flash mode, the cyclic and continuous modes which more closely simulate daylight exposure conditions. In each mode, the "fatigue resistance" parameter is computed. Moreover, the spectro-kinetic parameters of photomerocyanine species can be obtained from the kinetic mode. Due to the fully computer-controlled exptl. setup, the between-day reproducibility of the initial absorbance and fatigue resistance detns. are better than 4% and 6%, resp. Owing to the modular design of the

 $\mbox{\it menu-driven}$  software, written in C language, the app. is very flexible and easy to use.

- L35 ANSWER 51 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 119:84703 CA
- TI Development of an expert system for selection of experimental designs
- AU Olivero, Ramon A.; Seshadri, Sridhar; Deming, Stanley N.
- CS Lockheed Environmental Systems and Technologies Company, Las Vegas, NV, USA
- SO Analytica Chimica Acta (1993), 277(2), 441-53
- An expert system has been developed to assist chemists in the selection AΒ of exptl. designs for research projects. The system (named DXPERT) ranks thirteen types of exptl. designs according to their suitability for projects presented by users in an interactive session. categories included are factorial, response surface, sequential simplex optimization, simplex mixt., and statistical testing. A desirability index is assigned to each design alternative according to project characteristics (attributes). Characteristics are interpreted based on expert knowledge built into the system. DXPERT uses math. concepts to mimic features of human intuition and decision making. knowledge is represented by relevance factors (a no. between minus one and plus one) in a multiple-alternative multiple-attribute table. Relevance factors are interpreted as fuzzy values that represent the degree to which a design belongs to the set of suitable designs. formula for calcg. design desirabilities is based on fuzzy mathematics. For efficiency purposes, the order of the questions presented to the user is driven by a max. potential information gain algorithm. Design desirability indexes were found to be useful to researchers for the elimination of unsuitable designs and concn. of further efforts in the most applicable designs. A validation test was conducted with the participation of four other experts in the field.
- L35 ANSWER 57 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 118:80695 CA
- TI Application of an automated chemistry workstation to problems in **synthetic** chemistry
- AU Corkan, L. Andrew; Plouvier, Jean Christophe; Lindsey, Jonathan S.
- CS Dep. Chem., Carnegie Mellon Univ., Pittsburgh, PA, 15213, USA
- SO Chemometrics and Intelligent Laboratory Systems (1992), 17(1), 95-105
- AB An automated chem. workstation is applied to problems in the synthetic chem. of porphyrins. A factorial design study (16 expts., 96 data points) was performed to examine the role of catalyst and reactant concns. on porphyrin yield. Four expts. could be scheduled to run concurrently; all sixteen expts. were completed in less than 1 day of workstation time. The response surface from this expt. shows the conditions for achieving the highest yield. A simplex optimization was performed over the same reaction space, requiring fewer expts. to arrive at the optimal reaction parameters. A strategic search was performed to screen a list of reagents for catalytic activity. The effective concn. range of each catalyst was surveyed by systematic modification of an ongoing reaction. By terminating reactions when a

yield threshold was surpassed or when the entire concn. range had been spanned, compds. with catalytic activity and their effective concn. ranges were identified with minimal experimentation. Automated chem. workstations of this type should yield rapid accelerations in scientific research.

- L35 ANSWER 58 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 118:41343 CA
- TI Experiment manager software for an automated chemistry workstation, including a scheduler for parallel experimentation
- AU Corkan, L. Andrew; Lindsey, Jonathan S.
- CS Dep. Chem., Carnegie Mellon Univ., Pittsburgh, PA, 15213, USA
- SO Chemometrics and Intelligent Laboratory Systems (1992), 17(1), 47-74
- A review, with 10 refs., of an automation system with the ability to AΒ work relentlessly, precisely, strategically, and autonomously in pursuit of scientific goals. Some years' work has been aimed at developing the hardware and software architecture for an automated workstation. The workstation is designed for microscale experimentation in relatively clean domains of synthetic chem. the software system is presented. An expt. manager software package has been developed that provides for compn. of exptl. plans, controls all aspects of automated experimentation, and manages the data. manager software is comprised of open-loop and closed-loop expt. planners bundled together with supporting features for timing, scheduling, material data bases, resource management, automated startup procedures, running display, data handling, configuration options, and maintenance operations. The expt. planners provide a menu-driven user interface, editing features, and a modular set of procedural events with which diverse exptl. protocols can be composed. overarching objective has been to achieve versatility in expt. planning and still maintain access to the power that parallelism can confer in experimentation. The second theme is a description of approaches to performing expts. in parallel. Parallelism originates chiefly through simultaneous processing of samples at semi-autonomous hardware modules, at the user interface, and through the use of a scheduler. scheduler takes as input a set of exptl. protocols, establishes a sequence of the protocols, and interleaves the exptl. protocols without altering the relative times of the procedures within each protocol. Exptl. throughput can be increased by up to ten-fold by this approach.
- L35 ANSWER 60 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 117:233050 CA
- TI Microcomputers: bringing the process control laboratory into the classroom
- AU De Vaal, P. L.
- CS Dep. Chem. Eng., Univ. Pretoria, Pretoria, S. Afr.
- SO South African Journal of Chemical Engineering (1992), 4(1), 16-26
- AB A discussion is presented on the use of computers in the process control lab. in the Department of Chem. Engineering at the University of Pretoria. Examples are given of available software and of exptl. test app. and how these are used for educational and training purposes.

- ANSWER 65 OF 184 CA COPYRIGHT 2004 ACS on STN L35
- AN 115:138732
- Calo-Visco-Densimeter: a multipurpose apparatus TI
- Zhech, Helmut; Knapp, Helmut ÀU
- BASF, Ludwigshafen, D-6700, Germany CS
- Chemical Engineering & Technology (1991), 14(2), 109-13 SO
- A multipurpose app. was designed, built, and operated with the aim of AΒ the simultaneous detn. of enthalpy of mixing, heat capacity, kinematic viscosity, and d. of liq. mixts. at 240-370 K and 1-25 bars. operation of the app. was monitored by a personal computer. The design and the exptl. results were presented.
- ANSWER 84 OF 184 CA COPYRIGHT 2004 ACS on STN L35
- AN107:58093 CA
- The tin/lead solid/liquid phase diagram: a computer-controlled ΤI experiment
- Williams, Kathyrn R.; Eyler, John R.; Colgate, Samuel O. ΑU
- Univ. Florida, Gainesville, FL, 32611, USA CS
- Journal of Chemical Education (1987), 64(6), 499-500 SO
- A computer-interfaced detn. of the the solid-liq. phase diagram for the AB Sn/Pb system is described. Construction of the cooling curves (acquisition of thermocouple readings and subsequent plotting) can be performed with the aid of a microcomputer, leaving students adequate time to locate the break and arrest points, plot the phase diagram, and perform error anal.
- ANSWER 85 OF 184 CA COPYRIGHT 2004 ACS on STN L35
- ΑN 106:129098 CA
- Pulse radiolysis equipment: a setup for simultaneous multiwavelength TIkinetic spectroscopy
- Saran, Manfred; Vetter, Georg; Erben-Russ, Michael; Winter, Robert; ΑU Kruse, Alf; Michel, Christa; Bors, Wolf
- GSF Res. Cent., Inst. Strahlenbiol., Neuherberg, 8042, Fed. Rep. Ger. CS
- Review of Scientific Instruments (1987), 58(3), 363-8 SO
- A setup for pulse radiolysis expts. is described and consists of the AΒ following main components: an array of 15 photomultipliers attached to a spectrograph (allowing 4032 spectra to be recorded with a time resoln. of 500 ns/spectrum), self-regulating high-voltage supply for the photomultipliers, computer-controlled soln. mixing and dispensing system, and provisions for continuous dose variation of the FEBETRONaccelerator electron beam. The general performance of the system is discussed and construction or electronic details are given for special components.
- ANSWER 90 OF 184 CA COPYRIGHT 2004 ACS on STN L35
- AN 105:174871 CA
- Automated chemical synthesis. Part 4: Batch-type reactor automation TIand real-time software design
- Chodosh, Daniel F.; Kamholz, Kenneth; Levinson, Sidney H.; Rhinesmith, ΑU Robert

- CS Smith Kline and French Lab., Philadelphia, PA, 19101, USA
- SO Journal of Automatic Chemistry (1986), 8(3), 106-21
- AB The computer system, the temp. control system, the liq. trafficking system, the online chem. anal. system, the stirred-tank reactor, the software, and exptl. designs are discussed.
- L35 ANSWER 99 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 104:19150 CA
- TI Microcomputers in experimental materials science
- AU Gibeling, Jeffery C.
- CS Dep. Mech. Eng., Univ. California, Davis, CA, 95616, USA
- SO Comput. Usage Mater. Educ., Proc. Symp. (1985), Meeting Date 1984, 89-99. Editor(s): Liedl, Gerald L.; Sree Harsha, K. S. Publisher: Metall. Soc., Warrendale, Pa.
- The application of computer-based data acquisition and control techniques in senior undergraduate and 1st-year graduate level lab. courses is described. Areas of application include phase diagram detn., quant. metallog., Hall effect measurements, photocond. expts., creep testing, and conventional tensile testing.
- L35 ANSWER 101 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 102:51195 CA
- Automatically controlled measuring-computing complex for studying acoustic properties of matter in a wide range of parameters of state
- AU Otpushchennikov, N. F.; Tsyrenov, A. A.; Melikhov, Yu. F.; Mel'nikov, G. A.; Khanarin, V. S.
- CS USSR
- SO Ul'trazvuk Termodin. Svoistva Veshchestva (1983), 15-21. Editor(s): Otpushchennikov, N. F. Publisher: Kursk. Gos. Pedagog. Inst., Kursk, USSR.
- LA Russian
- AB An exptl. setup is presented for computer-controlled study of sound velocity and absorption in liqs., glasses, and polymers in wide temp., pressure, and frequency ranges.
- L35 ANSWER 114 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 96:93989 CA
- TI A semihierarchic **computer** network for data acquisition and **control** in staircase voltammetry
- AU Li, Chia-Yu; Barrett, Thomas H., Jr.; Lunney, David; Salt, Alger
- CS Dep. Chem., East Carolina Univ., Greenville, NC, 27834, USA
- SO Analytica Chimica Acta (1982), 134, 167-78
- AB A simple computer system designed for controlling voltammetric expts. was implemented. The system connects a Hewlett-Packard 2100A minicomputer with a Texas Instruments 9900 microcomputer by means of a bidirectional serial transmission link. Exptl. parameters are downloaded from the host to the satellite which supervises the expts. at a remote location. The data collected are transmitted back to the host at 2400 baud (data rate) for redn. and plotting. Through this division of labor, each system is used to its best advantage. Although the software developed is specifically for staircase voltammetry, the

system hardware is of general-purpose design which is suited for other types of pulse expt.

- L35 ANSWER 121 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 91:104620 CA
- TI Characterizing chemical systems with on-line computers and graphics
- AU Frazer, Jack W.; Rigdon, Lester P.; Brand, Hal R.; Pomernacki, Charles L.
- CS Lawrence Livermore Lab., Univ. California, Livermore, CA, 94550, USA
- SO Analytical Chemistry (1979), 51(11), 1739-47
- Incorporating computers and graphics on-line to chem. expts. and AB processes opens up new opportunities for the study and control of Systems having many variables can be characterized complex systems. even when the variable interactions are nonlinear, and the system cannot a priori be represented by numerical methods and models. sets of accurate data can be acquired rapidly, and then modeling and graphic techniques can be used to obtain partial interpretation plus design of further experimentation. The experimenter can thus comparatively quickly iterate between experimentation and modeling to obtain a final soln. A versatile computer-controlled app. was designed and characterized for chem. research, which incorporates on-line instrumentation and graphics. It can be used to det. the mechanism of enzyme-induced reactions or to optimize anal. methods. The app. can also be operated as a pilot plant to design control strategies. line graphics were used to display conventional plots used by biochemists and 3-dimensional response-surface plots.
- L35 ANSWER 135 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 82:158171 CA
- Instrumentation-oriented microcomputer. Extremely inexpensive data acquisition computer optimized for the automated laboratory
- AU Woodward, W. Stephen; Ridgway, Thomas H.; Reilley, Charles N.
- CS Kenan Lab. Chem., Univ. North Carolina, Chapel Hill, NC, USA
- SO Analyst (Cambridge, United Kingdom) (1974), 99(1185), 838-52
- AB The design and realization of a versatile, multiaccess, on-line digital computer system for instrumentation design, expt. control, and data interpretation were discussed.
- L35 ANSWER 152 OF 184 CA COPYRIGHT 2004 ACS on STN
- AN 73:83463 CA
- TI Interactive electronic analytical instrumentation based on computerized experimental design
- AU Jones, David Oldham; Perone, Sam P.
- CS Dep. of Chem., Purdue Univ., Lafayette, IN, USA
- SO Analytical Chemistry (1970), 42(11), 1151-7
- AB A hardware device has been constructed for use in stationary electrode polarographic expts. to automatically **optimize** the exptl. measurements. The basic instrumental parameters were detd. from expts. using a small lab. digital computer on-line. These parameters were then applied to the construction of a hardware device to perform the data anal. and collection, along with exptl. control operations. This allowed the

computing equipment to be freed for use in exptl. design while having a device to perform the routine analyses. The device, a simple hybrid computing system using medium-scale integrated circuits, was designed to interact with the expt. in realtime to modify the exptl. parameters and to provide a sample-oriented anal. The crit. exptl. design parameters elucidated by the software-oriented study are considered and the translation to hardware-oriented instrumentation is described. addn., the ease with which efficient lab. instrumentation can be designed and constructed as a result of the new integrated circuit electronics technology is demonstrated. The hardware device was evaluated exptl. in the lab. and the capabilities were compared to the software-generated experimentation. Results indicate that essentially identical anal. performance is achieved at greatly reduced complexity of operation and expense. In addn., the hardware-oriented instrumentation provides the capability for greater measurement

L35 ANSWER 175 OF 184 CA COPYRIGHT 2004 ACS on STN

AN 51:69432 CA

OREF 51:12566d-f

TI Apparatus for differential thermal analysis of lubricating greases

AU Cox, David B.; McGlynn, James F.

CS Socony Mobil Oil Co., Inc., Brooklyn, NY

SO Anal. Chem. (1957), 29, 960-3

An app. is described for the investigation of Li soap-oil systems by using differential thermal analysis. The app. operates automatically in controlling heating rates and recording all pertinent data and differential temps. Five different detns. may be run simultaneously. The design of the cell and thermopile permits high sensitivity and accuracy in locating phase changes. A differential thermal analysis of Li stearate revealed a previously unreported mesomorphic phase occurring between 225°and 229°. A more accurate phase diagram than previously reported was obtained of the Li stearate-n-hexadecane system.

L35 ANSWER 180 OF 184 CA COPYRIGHT 2004 ACS on STN

AN 41:11186 CA

OREF 41:2283f-i,2284a

TI Recording viscometer for starches

AU Kesler, C. C.; Bechtel, W. G.

CS Penick & Ford, Ltd., Cedar Rapids, Iowa

SO Journal of Industrial and Engineering Chemistry (Washington, D. C.) (1947), 39, 16-21

AB A continuous-recording viscometer for routine and research testing of starch products is described. Variations in cooking procedures which cause errors in  $\eta$  detns. are prevented through **automatic control** of **heating**, max. **temp**., rate of stirring, and loss of water by evapn.  $\eta$  is measured as the force which the paste exerts against a propeller driven through it at const. speed. A gear differential transmits the force to a dynamometer attached to the pen arm of the recorder. Interchangeable wts. on the dynamometer permit measurements in several ranges with equal sensitivity and without requiring recalibration of the

viscometer. Viscosities of one poise or more can be measured throughout any desired length of time, and during both heating and cooling periods. Variations between samples and differences between com. types of starch can be readily observed. Starches from different sources and of different kinds and degrees of modification give characteristically different curves which are of value in their identification and study. Results of 75 detns. in duplicate made by the same operator show the precision to be expected is within  $\square$ %. Where a single starch was tested by different operators in different labs. the degree of precision was approx.  $\square$ %. Studies of sensitivity with a concn. of 5% of unmodified corn starch as the control showed that to avoid errors due to measurement of samples the starch should be weighed to within  $\Box$ .05 g. and the water measured to within  $\pm 1$  ml. vol. of the starch soln. may be varied by as much as 100 ml. without affecting the readings. Design of the agitator permits measurement of  $\eta$  up to about 55 poises; this covers the range from 3% unmodified tapioca or potato starch or 4% unmodified corn starch to 6.5% tapioca or 9 to 10% cornstarch.

=> log y STN INTERNATIONAL LOGOFF AT 13:41:51 ON 19 OCT 2004